Medical

Scientific & Technical Report

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Referenze

1. Trautmann M, Royer H, Helm E, May W, Haller M. *Pseudomonas aeruginosa*: new insights into transmission pathways between hospital water and patients. *Filtration* 2004; 1(Supplement): 63-70
2. Blanc DS, Nahimana I, Petignant I, Wenger A, Bille J, Francioli P. Faucets as a reservoir of edemic *Pseudomonas aeruginosa* colonization/infections in intensive care units. *Intensive Care Med* 2004; 30: 1964-8
3. Warris A, Ustad P, Meis JFGM, Voss A, Verweij PE, Abrahamsen TG. Recovery of filamentous fungi from water in paediatric bone marrow transplant unit. *Journal of Hospital Infection* 2001; 47: 143-8
4. Sabria M, Yu VL. Hospital-acquired legionellosi: solutions for preventable infection. *Lancet Infect Dis.* 2002 2(6): 368-73
5. Patterson WJ, Hay J, Seal DV, McLuckie JC. Colonization of transplant unit water supplies with *Legionella* and protozoa: precautions required to reduce the risk of legionellosis. *Journal of Hospital Infection* 1997; 26: 417-425
6. Darelid j, Bergtsson L, Gastrin B et al. An outbreak of Legionnaires Disease in a Swedish Hospital. *Scandinavian J. Infectious Disease.* 1994; 26: 417-425
7. Marrie TJ, Haldne D, Bezanson G, Peppard R. Each water outlet is a unique ecological niche for *Legionella pneumophila*. *Epidemiology Infection* 1992; 108: 261-270
8. Goetz Am, Stout JE, Jacobs SL et al. Nosocomial legionnaires disease discovered in community hospitals following cultures of the water system: seek and ye shall find. *American j Infection Control* 1998; 26: 8-11
9. From the centers for Disease Control Prevention. Transmission of nosocomial Legionnaires disease. *J American Medical Assoc* 1997; 277: 1927-8
10. Ricci ML et al. A 4-year experience of ClO₂-treatment for the control of *Legionella* contamination in a hospital water system. 20th Meeting EWGLI, Rome, P18 16/17th May 2005
11. Hosein IK, Hill DW, Tan TY et al. Point-of-care controls for nosocomial legionellosis combined with chlorine dioxide potable water decontamination: a two-year survey at a Weksh teaching hospital. *Journal of Hospital Infection* 2005; 61: 100-106
12. Exener M, Kramer A, Lajoie L, gebel J, Engelhart S, Hartemann P. Prevention and control of health care-associated waterborne infections in health care facilities. *Am J infect Control* 2005; 33: S26-40

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I part numbers descritti sono protetti da copyright di Pall Europe Limited.
Pall e Pall-Aquasafe sono Marchi di Fabbrica di Pall Corporation.
Filtration. Separation. Solution. è Marchio di servizio di Pall Corporation.
© 2008 Pall Europe Limited. PMED/1M/DBD/1.2006