

Pall Medical
A division of Pall Europe Limited
Walton Road
Farlington
Portsmouth PO6 1TD, England

+44 (0)23 9230 3303 phone
+44 (0)23 9230 3541 main office fax
+44 (0)23 9230 3324 sales office fax
www.pall.com

14th July 2009

To whom it may concern

RE: *Pall-Aquasafe* Water Filter Field Evaluation Reports


Pall-Aquasafe Water Filters have been validated using liquid microbial challenge tests (*Brevundimonas diminuta* (ATCC 19146)) using an industry standard test method for 0.2 µm sterilising grade filters to $\geq 10^7$ CFU/cm² of effective filtration area.

In addition Pall often completes Field Evaluation Reports which contain data generated within health care facilities.

Although Pall provide guidance and filter usage protocols, the selected health care facilities are independent to Pall and microbiological test samples generated from field evaluations are sent to accredited independent laboratories for analysis.

Examples of such Field Evaluation Reports can be requested via your local Pall Sales Office.

Yours sincerely,



Vanessa Merefield
Life Sciences Scientific & Laboratory Services Manager

Pall Italia s.r.l.
Via Emilia, 26
20090 Buccinasco - MI - Italia
Telefono +39 02 488870 1
Fax +39 02 4880014
www.pall.com

14 luglio 2009

Oggetto: Report di Valutazione in campo dei Filtri per acqua Pall-Aquasafe


Filtri per acqua Pall-Aquasafe sono stati validati utilizzando il test di “challenge” microbico liquido” (*Brevundimonas diminuta* (ATCC 19146)), seguendo la metodica standard dell’industria per i filtri di grado sterilizzante da 0.2µm, alla concentrazione $\geq 10^7$ UFC/cm² di superficie effettiva di filtrazione.

A completamento, Pall spesso aggiunge nei Report Tecnici delle Valutazioni in Campo tutti i dati generati all’interno delle strutture sanitarie.

Pall fornisce dunque guide e procedure per l’utilizzo dei filtri; le selezionate strutture sanitarie sono indipendenti da Pall e i test microbiologici dei campioni generati durante queste “valutazioni in campo” vengono inviati a indipendenti laboratori accreditati.

Esempi di Report Tecnici delle Valutazioni in Campo possono essere richiesti agli uffici Pall.

Sentitamente,


Dott.ssa Cinzia Quarti
Responsabile Tecnico e Marketing
Pall Italia S.r.l.
Divisione Lifesciences